

Figure 1

1/10

1 FPTIPLSRFL DNASLRAHRL HQLAFDITYQE FEEAYIPKEQ KYSFLQNPQT  
51 SLCFSESIPT PSNREETQK SNLELLRISL LLIQSWLEPV QFLRSVFANS  
101 LVGASDSNV YDLLK DLEEGIQ TLMGRLED GSPRTGQIFK QTYSKFDTNS  
151 HNDDALLKNY GLLYCFRKDM DKVETFLRIV QCRSVEGSCG F

Figure 2

2/10

TTCCCAACCATTCCTTATCCAGGCTTTTGGACAACGCTATGCTCCGCGCCCATCGTCTGCACCAGCTGGC  
CTTTGACACCTACCAGGAGTTTGAAGAAGCCTATATCCCAAAGGAACAGAAATTCATTCTGCAGAAC  
CCCCAGACCTCCCTCTGTTTCTCAGAGTCTATTCCGACACCCTCCAACAGGGAGGAAACACAACAGAAAT  
CCAACCTAGAGCTGCTCCGCATCTCCCTGCTGCTCATCCAGTCGTGGCTGGAGCCCGTGCAGTTCCTCAGG  
AGTGTCTTCGCCAACAGCCTGGTGTACGGCGCCTCTGACAGCAACGTCTATGACCTCCTAAAGGACCTAG  
AGGAAGGCATCCAAACGCTGATGGGGAGGCTGGAAGATGGCAGCCCCGGACTGGGCAGATCTTCAAGC  
AGACCTACAGCAAGTTCGACACAACTCACACAACGATGACGCACTACTCAAGAACTACGGGCTGCTCTA  
CTGCTTCAGGAAGGACATGGACAAGGTCGAGACATTCCTGCGCATCGTGCAAGTGCCGCTCTGTGGAGGGC  
AGCTGTGGCTTC


Figure 3

3/10

MDLWQLLLTLALAGSSDAFSGSEATAAILS RAPWSLOS VNPGLKTNSSKEPKFTKCRSPER  
ETF SCHWTDEVHHGTKNLGPIOLFYTRRNTQEWTOEWKECPDYVSAGENSCYFNSSFTSI  
WIPYCIKLT SNGGT VDEKCF SVDEIVQDPPIALNWTL LNVS LTGIHADIQVRWEAPRNADI  
OKGWMVLEYELO YKEVNETKWKMMDPILTTSVPVYSLKVDKEYEVRVR SKORNSGNYG  
EFSEVLYVTLPOMS QFTCEEDFYFPWLLIIIFGIFGLTVMLFVFLFSKQQRIKMLILPPVPVPK  
IKGIDPDLLKEGKLEEVNTILAIHDSYKPEFHSDDSWVEFIELDIDEPDEKTEESD TDRLSSD  
HEKSHSNLGVKDGDSGRTSCCEPDILETDFNANDIHEGTSEVAQPQRLKGEADLLCLDQKN  
QNNSPYHDACPATQQPSV IQAEKNKPQPLPTEGAESTHQA AHIQLSNPSSLSNIDFYAQVSD  
ITPAGSVVLSPGQKNKAGMSQCDMHPEMVSLCQENFLMDNAYFCEADAKK CIPVAPHIKV  
ESHIQPSLNQEDIYITTESLT TAAGRPGTGEHVPGSEMPVPDYTSIHIVQSPQGLILNATALPL  
PDKEFLSSCGYVSTDQLNKIMP

Figure 4  
Strategy

4/10



- 1) PCR1 with forward (GHPermN0x+) and reverse (GHPermL01-) primers.
- 2) Purify PCR product.
- 3) PCR2 with product from PCR1 and reverse primer (GHPermC0x-).
- 4) Purify PCR product.
- 5) Digest PCR product and suitable vector using *Bam*HI and *Not*I.
- 6) Ligate PCR product into vector.


**Primer Pairs/PCRs**

	First PCR		Second PCR	
	Forward	Reverse	Forward	Reverse
GHCP02	GHPermN01+	GHPermL01-	PCR1 product	GHPermC02-
GHCP03	GHPermN01+	GHPermL01-	PCR1 product	GHPermC03-
GHCP04	GHPermN04+	GHPermL01-	PCR1 product	GHPermC01-
GHCP05	GHPermN04+	GHPermL01-	PCR1 product	GHPermC02-
GHCP06	GHPermN06+	GHPermL01-	PCR1 product	GHPermC06-
GHCP07	GHPermN07+	GHPermL01-	PCR1 product	GHPermC07-

BEST AVAILABLE COPY

Figure 5

5/10



Growth Hormone

1 TTCCAACCATTCCTTATCCAGGCTTTTGGACAACGCTAGTCTCCGCGC  
 51 CCATCGTCTGCACCAGCTGGCCTTTGACACCTACCAGGAGTTTGAAGAAG  
 101 CCTATATCCCAAAGGAACAGAAGTATTCAATTCCTGCAGAACCCCCAGACC  
 151 TCCCTCTGTTTCTCAGAGTCTATTCCGACACCCTCCAACAGGGAGGAAAC  
 201 ACAACAGAAATCCAACCTAGAGCTGCTCCGCATCTCCCTGCTGCTCATCC  
 251 AGTCGTGGCTGGAGCCCGTGCAGTTCCCTCAGGAGTGTCTTCGCCAACAGC  
 301 CTGGTGACGGCGCCTCTGACAGCAACGTCTATGACCTCCTAAAGGACCT  
 351 AGAGGAAGGCATCCAAACGCTGATGGGGAGGCTGGAAGATGGCAGCCCCC  
 401 GGACTGGGCAGATCTTCAAGCAGACCTACAGCAAGTTCGACACAACTCA  
 451 CACAACGATGACGCACTACTCAAGAACTACGGGCTGCTCTACTGCTTCAG  
 501 GAAGGACATGGACAAGGTCGAGACATTCTGCGCATCGTGCAGTGCCGCT  
 551 CTGTGGAGGGCAGCTGTGGCTTC


GH\_CP01

1 ATCCAAACGCTGATGGGGAGGCTGGAAGATGGCAGCCCCCGGACTGGGCA  
 51 GATCTTCAAGCAGACCTACAGCAAGTTCGACACAACTCACACAACGATG  
 101 ACGCACTACTCAAGAACTACGGGCTGCTCTACTGCTTCAGGAAGGACATG  
 151 GACAAGGTCGAGACATTCTGCGCATCGTGCAGTGCCGCTCTGTGGAGGG  
 201 CAGCACCATTCCCTTATCCAGGCTTTTGGACAACGCTAGTCTCCGCGCCC  
 251 ATCGTCTGCACCAGCTGGCCTTTGACACCTACCAGGAGTTTGAAGAAGCC  
 301 TATATCCCAAAGGAACAGAAGTATTCAATTCCTGCAGAACCCCCAGACCTC  
 351 CCTCTGTTTCTCAGAGTCTATTCCGACACCCTCCAACAGGGAGGAAACAC  
 401 AACAGAAATCCAACCTAGAGCTGCTCCGCATCTCCCTGCTGCTCATCCAG  
 451 TCGTGGCTGGAGCCCGTGCAGTTCCCTCAGGAGTGTCTTCGCCAACAGCCT  
 501 GGTGTACGGCGCCTCTGACAGCAACGTCTATGACCTCCTAAAGGACCTAG  
 551 AG

BEST AVAILABLE COPY

Figure 6

6/10



BEST AVAILABLE COPY

Figure 7

7/10

GH121 to GH118






Figure 8

8/10



BEST AVAILABLE COPY



TABLE 1

9/10

Primer Name	Sequence (5' → 3')	Length
GH_CP01For	gctaggatccaacccttATCCAAACGCTGATGG	33
GH_CP01Link	tggataaggggaatggtGCTGCCCTCCACAGAG	32
GH_CP01Rev	gtcaactggtcagcggccgccCTCTAGGTCCTTTAGGAG	39

TABLE 2

10/10

<b>Primers (all 5'→3')</b>	
<b>Old Termini Linkers</b>	
GHPermL01+	ctctgtggagggcagcACCATTCCCTTATCCA (32)
GHPermL01-	tggataaggggaatggtGCTGCCCTCCACAGAG (32)
<b>GHCP01 (121→118)</b>	
GHPermN01+	gctaggatccaacccttATCCAAACGCTGATGG (33)
GHPermC01- (39)	gtcaactggtcagcggccgcccCTCTAGGTCCTTTAGGAG
<b>GHCP02 (121→119)</b>	
GHPermN01+	gctaggatccaacccttATCCAAACGCTGATGG (33)
GHPermC02- (40)	gtcaactggtcagcggccgcccTCCTCTAGGTCCTTTAGG
<b>GHCP03 (121→120)</b>	
GHPermN01+	gctaggatccaacccttATCCAAACGCTGATGG (33)
GHPermC03-	gtcaactggtcagcggccgcccGCCTTCCTCTAGGTCC (37)
<b>GHCP04 (120→118)</b>	
GHPermN04+	gctaggatccaacccttGGCATCCAAACGCTGATGG (36)
GHPermC01- (39)	gtcaactggtcagcggccgcccCTCTAGGTCCTTTAGGAG
<b>GHCP05 (120→119)</b>	
GHPermN04+	gctaggatccaacccttGGCATCCAAACGCTGATGG (36)
GHPermC02- (40)	gtcaactggtcagcggccgcccTCCTCTAGGTCCTTTAGG
<b>GHCP06 (102→100)</b>	
GHPermN06+ (38)	gctaggatccaacccttGTGTACGGCGCCTCTGACAGC
GHPermC06- (41)	gtcaactggtcagcggccgcccGCTGTTGGCGAAGACACTCC
<b>GHCP07 (132→130)</b>	
GHPermN07+	gctaggatccaacccttAGCCCCCGGACTGGGCAG (35)
GHPermC07- (40)	gtcaactggtcagcggccgcccATCTTCCAGCCTCCCCATC